Blood & Urine Kit Instruction Sheet

(Combination Blood & Urine Specimen Collection Kit)

For Use by Medically Trained Personnel Only

IMPORTANT! PLEASE READ INSTRUCTIONS CAREFULLY BEFORE PROCESSING

Contents:

One:	2-part identification (ID) form with bar code labels and tamper-evident tape	One:	Alcohol prep
One:	Prep App Info	One:	Sterile gauze pad
Two:	Non-Latex gloves	One:	Bandage
One:	EDTA tube (purple top)	One:	Urine collection cup with temperature strip affixed
One:	Tiger Top (tube)	One:	Blue urine tube & cap
One:	Serum transfer tube (clear tube and cap)	One:	Green urine tube & cap with preservative tablet
One:	Blood Collection Needle	One:	Leak-resistant plastic bag
One:	Needle holder	One:	Absorbent sheet
		One:	Approved shipping box

Notice:

INTENDED USE: This kit is collected blood and/or urine samples for diagnostic laboratory testing.

WARNING – USE UNIVERSAL PRECAUTIONS: All human blood and urine samples should be handled and processed as if potentially infectious.

WASTE DISPOSAL INSTRUCTIONS:

1. Dispose of all needles in an approved sharps container.
2. Dispose of medical waste (gloves, wipes, etc.) by approved methods for medical waste.

PREPARING THE LAB ID / CONSENT FORM & TAMPER-EVIDENT TAPE

1. Complete all information requested on the ID form except:
 - Urine temperature
 - Examiner chain-of-custody sections
 - Proposed insured chain of custody sections
2. Have proposed insured sign on the tamper-evident tape.
3. Affix barcode labels with the proposed insured's full name to all the tubes, except the blue tube. The BLUE urine tube will be labeled with the bar coded tamper-evident tape after collection.

INSTRUCTIONS FOR COLLECTION OF URINE SPECIMEN

1. Examiner should wear gloves.
2. Instruct the proposed insured to void into the collection cup and return it to you immediately (cup with temperature strip affixed). Urine level should be above the temperature strip for accurate reading.
3. Read temperature strip on collection cup.
 - Wait two (2) minutes, then read the temperature (brightest portion of the line) and record this number in the "URINE TEMP" space on the ID form.
 - If the temperature strip does not activate with colors, record temperature as "OTR" (outside the range).
4. Pour urine from collection cup into the two (2) urine tubes.
 - Fill each tube approximately ¾ full. Do not overflow.
 - Cap the tubes securely.
 - Center the signed tamper-evident tape over the cap of the BLUE urine tube.
 - Place both tubes in plastic bag.
 - Dispose of the urine collection cup by approved methods for medical waste.

INSTRUCTIONS FOR COLLECTION OF BLOOD SPECIMEN

1. Draw blood from the proposed insured into the evacuated tube.
2. Immediately after the draw, gently invert the EDTA (purple-top) tube 10 times to mix with anti-coagulant, and return the tube to the kit.
3. Immediately after the draw, gently invert serum/gel (tiger top) tube 6-8 times allowing blood to flow completely from one end to the other. Do not shake: vigorous mixing may cause hemolysis. Let the serum/gel (tiger top) stand upright until the blood has clotted (at least 15 minutes, not longer than 30 minutes). Do not lay tube flat. After clotting, the serum/gel (tiger top) tube must be centrifuged for 15 minutes at 3200 RPMs. Inspect serum/gel tubes after centrifugation to ensure that a complete barrier has formed between coagulum and serum. Transfer serum to the serum-transfer (CLEAR) tube and cap. Serum must be red-cell free. Place centrifuged serum/gel (tiger top) tube in the plastic bag.
4. Follow packaging instructions for express shipping.

IF BLOOD IS DRAWN ON MOBILE BASIS WITHOUT IMMEDIATE ACCESS TO CENTRIFUGE:

1. Keep the serum/gel (tiger top) tube upright and store in cooler.
2. Centrifuge within 4 hours.
3. Transfer serum to the serum transfer (CLEAR) tube and cap.
4. DO NOT FREEZE.

AFTER SPECIMEN COLLECTION

1. Complete the chain-of-custody sections for examiner and proposed insured.
2. Provide "Important Applicant Brochure" to the proposed insured.

PACKAGING & SHIPPING INSTRUCTIONS

1. Place tubes and vials upright in leak-resistant bag.
2. Do not remove absorbent sheet from plastic bag.
3. Apply firm pressure over orange zipper fastener.
4. Fold ID around the plastic bag and place in box. Fold all additional paperwork and place adjacent to the plastic bag in the box.
5. Close box.
6. Affix the designated bar code located on the label pod in the space provided on the outside of the box.
- 7a. Insert box in Lab pack.
- 7b. Fill out airbill address to:
 - ExamOne
 - 10101 Renner Blvd.
 - Lenexa, KS 66219, U.S.A.
- 7c. Mark the airbill: Bill to: Receiver

CAUTION

The urine preservative tablet in this container is NOT meant for human consumption. The tablet is to remain in the container with the urine.

Ordinarily, the level of preservative is low enough that harm does not follow ingestion. However, if ingested, we strongly recommend that you immediately consult a poison control center, your family physician or the nearest hospital emergency room.

The Cargille Stabilur Tablet in the green tube contains the following hazardous ingredients.

Mercuric Oxide	Hexamethylene Tetramine
Benzoic Acid	Bicarbonate of Soda
Mono Potassium Phosphate	Sodium Benzoate

BD VACUTAINER™ ECLIPSE™ BLOOD COLLECTION NEEDLE

- 1) Holding both colored shields, twist and remove white shield
- 2) Screw on holder
- 3) Rotate safety shield back
- 4) Twist and pull needle shield straight off
- 5) Perform venipuncture
- 6) Firmly push forward on the safety shield, lock into place and inspect

Things to Remember:

- Handle all biologic samples and blood collection sharps, lancets, needles, luer adapters and blood collection sets according to the policies and procedures of your facility.
- Obtain appropriate medical attention in the event of any exposure to biologic samples (for example, through a puncture injury) since they may transmit viral hepatitis, HIV (AIDS), or other infectious diseases.
- Utilize any built-in needle protectors in the blood collection device.
- Discard needle in the nearest sharps collector. Follow your institution's policy for safe disposal of all medical waste.

DO NOT USE ANY USED OR UNUSED BLOOD COLLECTION NEEDLES TO EXAMONE.

If you have kit concerns or order problems, please call our Kit Department at 1-800-952-2350.

Quest, Quest Diagnostics, ExamOne, any associated logos, and all associated Quest Diagnostics registered or unregistered trademarks are the property of Quest Diagnostics. All third-party marks - © and ™ - are the property of their respective owners. © 2015 Quest Diagnostics Incorporated. All rights reserved.

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